


# THORACIC Instruments




Bone and Rib Shears, Sternum Chisel



GLUCK  
28-002-19  
19 cm / 7 1/2"




28-003-22  
22 cm / 8 3/4"




28-007-01  
For the Right Hand  
28-007-02  
For the Left Hand

BRUNNER  
28-007-01 – 28-007-02  
32 cm / 12 1/2"




GIERTZ  
28-014-25  
25 cm / 10"





LEBSCHKE  
28-019-90  
25 cm / 10"

Rib Spreaders



HAIGHT


	a	b	c
28-023-09	3 cm	3 cm	9 cm / 3 1/2"



TUFFIER


	a	b	c
28-022-16	4.5 cm	4.5 cm	16 cm / 6 1/4"

Rib Spreaders




FINOCHIETTO – BABY

	a	b	c
<b>28-021-75</b>	15 mm	20 mm	75 mm / 2"



FINOCHIETTO – HAIGHT - BABY

	a	b	c
<b>28-027-10</b>	31 mm	31 mm	100 mm / 4"




1/2

FINOCHIETTO

	a	b	c
<b>28-032-01</b>	30 mm	40 mm	150 mm / 6"
<b>28-032-02</b>	45 mm	65 mm	200 mm / 8"
<b>28-032-03</b>	75 mm	65 mm	250 mm / 10"


Rib Spreaders



BURFORD-FINOCHIETTO  
 with 2 pairs of blades

	a	b	c
28-033-20	65 mm 47 mm	62 mm 62 mm	200 mm / 8"
28-033-25	65 mm 47 mm	62 mm 62 mm	250 mm / 10"
28-033-30	65 mm 47 mm	62 mm 62 mm	300 mm / 12"

1/2



FINOCHIETTO-BABY

	a	b	c
28-020-90	18 mm	15 mm	90 mm / 3 1/2"

1/2


FINOCHIETTO

	a	b	c
28-034-01	30 mm	45 mm	150 mm / 6"
28-034-02	46 mm	62 mm	190 mm / 7 1/2"
28-034-03	80 mm	65 mm	190 mm / 7 1/2"


Rib Contractors, Raspatories, Lung Spatula







BAILEY  
 28-041-16  
 16 cm / 6 1/4"




BAILEY-GIBBON  
 28-042-16  
 16 cm / 6 1/4"




DOYEN  
 28-045-01 – 28-046-02  
 17 cm / 6 1/4"




ALLISON  
 28-050-27  
 27 cm / 10 1/4"  
 28-050-32  
 32 cm / 12 1/2"



Lung Grasping Forceps



LOVELACE  
28-056-19  
19 cm / 7 1/2"



LOVELACE  
28-057-19  
19 cm / 7 1/2"



Bronchus Forceps

